## Single Package (Front side)

## Single Package (Backside)



## **BACKYAK99** Refreshing MASK **Product Features** Functional nose-supporting system **BACKYAK99 SALT FILTER** Excellent fit due to ergonomic design Comfortable breathing Wearing Method Unfold the mask Precautions for Use Do not use in enclosed areas with oxygen concentrations below 18%. Do not wear it over wrapping the respiratory system with towels, tissues, etc. Do not wash and use. (filtering effect gets weak). Do not crush or deform the face of mask. Hyou are uncomfortable with breathing due to wearing masks, such as pregnances. women, respiratory and cardiovascular patients, children and the elde stop using them and consult with experts such as doctors if necessary **Product Notation** Product Name acrocare BACKYAK99 refreshing Mask This product may be exchanged or refunded in a resolution standards announced by the Fair Trad MADE IN KOREA